



Product Datasheet

HNF4A antibody (orb329675)


Description	Rabbit polyclonal antibody to HNF4A
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Porcine, Rat, Sheep, Zebrafish
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human HNF4A
Target	HNF4A
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: RGQAATPETPQSPGGSGSEPYKLLPGAVATIVKPLSAIPQPTITKQEV
Clonality	Polyclonal
MW	50kDa
Uniprot ID	P41235
NCBI	NP_787110
Expiration Date	12 months from date of receipt.



Western blot analysis of OVCAR-3 cell ly...



Immunohistochemical staining of human Sm...



Host: Rabbit, Target: HNF4A, Positive co...