


Product Datasheet

FOXP2 antibody (orb329741)

Description	Rabbit polyclonal antibody to FOXP2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human FOXP2
Target	FOXP2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: SGLKSPKSSDKQRPLQVPVSVAMMTPQVITPQQMQQILQQVLSPOQLQA
Clonality	Polyclonal
MW	80kDa
Uniprot ID	O15409
NCBI	NP_055306
Expiration Date	12 months from date of receipt.



Rabbit Anti-FOXP2 Antibody, Paraffin Emb...



WB Suggested Antibody Titration: 0.2-1 u...



WB Suggested Anti-FOXP2 Antibody, Titrat...