


Product Datasheet

VDAC2 antibody (orb329811)


Description	Rabbit polyclonal antibody to VDAC2
Species/Host	Rabbit
Reactivity	Canine, Equine, Goat, Guinea pig, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human VDAC2
Target	VDAC2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: VKLDVKTKSCSGVEFSTSGSSNTDTGKVTGTLETKYKWCEYGLTFTEKWN
Clonality	Polyclonal
MW	31kDa
Uniprot ID	P45880
NCBI	NP_003366
Expiration Date	12 months from date of receipt.



Western blot analysis of human Fetal Bra...



Western blot analysis of 293T cell lysat...



Host: Rabbit,
Target Name:
VDAC2,
Sample...