



Product Datasheet

Ppargc1a antibody (orb330094)

Description	Rabbit polyclonal antibody to Ppargc1a
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of mouse Ppargc1a
Target	Ppargc1a
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: QEIRAE LNKHFGHPCQAVFDDKSDKTS E LRDGDFSNEQFSKLPVFINSGL
Clonality	Polyclonal
MW	90kDa
Uniprot ID	O70343
NCBI	NP_032930
Expiration Date	12 months from date of receipt.



Western blot analysis of SP2/0 cell lysa...



Western blot analysis of SP2/0 tissue us...



Western blot analysis of Transfected Mel...