



Product Datasheet

SLC5A7 antibody (orb330358)

Description	Rabbit polyclonal antibody to SLC5A7
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C-terminal region of human SLC5A7
Target	SLC5A7
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: ILVKNENIKLDELALVKPRQSM T L S S T F T N K E A F L D V D S S P E G S G T E D N L
Clonality	Polyclonal
MW	63 kDa
Uniprot ID	Q9GZV3
NCBI	NP_001291934.1
Expiration Date	12 months from date of receipt.



25 ug of the indicated Human whole cell ...



Host: Rabbit,
Target Name:
SLC5A7,
Sampl...