



Product Datasheet

Nlgn3 antibody (orb330582)


Description	Rabbit polyclonal antibody to Nlgn3
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N-terminal region of Rat Nlgn3
Target	Nlgn3
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: KQGEDLADNDGDEDEDIRDSGAKPVMVYIHGGSYMEGTGNMIDGSVLASY
Clonality	Polyclonal
MW	93kDa
Uniprot ID	Q62889
NCBI	NP_599163
Expiration Date	12 months from date of receipt.



Western blot analysis of rat Testis tiss...



Host: Rabbit,
Target Name:
NLGN3,
Sample...



Host: Rabbit,
Target Name:
Nlgn3, Sample...