



Product Datasheet

HYAL1 antibody (orb330624)


Description	Rabbit polyclonal antibody to HYAL1
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human HYAL1
Target	HYAL1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: WNANTQWCLERHGVDVDVSVFDVVANPGQTFRGPDMTIFYSSQLGTYPY
Clonality	Polyclonal
MW	45kDa
Uniprot ID	Q12794
NCBI	NP_695014
Expiration Date	12 months from date of receipt.



Western blot analysis of human HepG2 tis...



Immunohistochemical staining of human Br...



Immunohistochemical staining of human Li...