



Product Datasheet

UXT antibody (orb330628)

Description	Rabbit polyclonal antibody to UXT
Species/Host	Rabbit
Reactivity	Canine, Equine, Goat, Guinea pig, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human UXT
Target	UXT
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: MATPPKRRAVEATGEKVLRYETFISDVLQRDLRKVLDHRDKVYEQLAKYL
Clonality	Polyclonal
MW	18 kDa
Uniprot ID	Q9UBK9
NCBI	NP_705582
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cell lysa...



25 ug of the indicated Human whole cell ...



WB Suggested Anti-UXT Antibody Titration...