



Product Datasheet

HK2 antibody (orb330743)


Description	Rabbit polyclonal antibody to HK2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human HK2
Target	HK2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: GTEHGEFLALDLGGTNFRVLWVKVTDNGLQKVEMENQIYAIPEDIMRGSG
Clonality	Polyclonal
MW	102kDa
Uniprot ID	P52789
NCBI	NP_000180
Expiration Date	12 months from date of receipt.



Host: Rabbit, Target:
HK2, Positive cont...



Immunohistochemistry
with Human Skeletal...



Lanes: 1: 50 ug HEP3B
lysate, 2: 50 ug H...