



Product Datasheet

GRN antibody (orb330767)


Description	Rabbit polyclonal antibody to GRN
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of Human GRN
Target	GRN
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: QAVCCEDHIHCCPAGFTCDTQKGTCEQGPVWMEKAPAHLSLPDPQAL
Clonality	Polyclonal
MW	64 kDa
Uniprot ID	P28799
Expiration Date	12 months from date of receipt.



25 ug of the indicated Human whole cell ...



50 ug of HEK293T WT and GRN KO lysates w...



Host: Rabbit, Target Name: GRN, Sample T...