


Product Datasheet

SERPINI2 antibody (orb330880)

Description	Rabbit polyclonal antibody to SERPINI2
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human SERPINI2
Target	SERPINI2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: LPRFKVEQKVDKFDVLYSLNITEIFSGGCDLSGITDSSEVYVSQVTQKVF
Clonality	Polyclonal
MW	46kDa
Uniprot ID	O75830
NCBI	NP_006208
Expiration Date	12 months from date of receipt.



Western blot analysis of COLO205 cell ly...



WB Suggested Anti-SERPINI2 Antibody Titr...