

Product Datasheet

Ip6k1 antibody (orb331181)

Description	Rabbit polyclonal antibody to Ip6k1
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	Ip6k1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: YPYVESETVEQDDTPEREQPRRKHSRRSLHRSGSGSDHKEEKASLSFETS
Clonality	Polyclonal
MW	49kDa
Uniprot ID	Q6PD10
NCBI	NP_038813
Expiration Date	12 months from date of receipt.



WB Suggested
Anti-*Ip6k1*
Antibody,
Titrat...