



Product Datasheet

MTNR1B antibody (orb331264)

Description	Rabbit polyclonal antibody to MTNR1B
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	MTNR1B
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: SFANCCEAGGWAVRPGWWSGAGSARPSRTPRPPWVAPALSAVLIVTTAVDV
Clonality	Polyclonal
MW	40 kDa
Uniprot ID	P49286
NCBI	NP_005950
Expiration Date	12 months from date of receipt.



25 ug of the indicated Human whole cell ...



WB Suggested Anti-MTNR1B Antibody, Titra...