

Product Datasheet

Connexin 37 antibody (orb33239)

Description Rabbit polyclonal antibody to Connexin 37 belongs to the connexin family that is expressed in multip...

Species/Host Rabbit

Reactivity Mouse

Conjugation Unconjugated

Tested Applications IHC-P, WB

Immunogen Synthetic Peptide

Target GJA4

Form/Appearance Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal



MW 37414

Uniprot ID [P35212](#)

NCBI [NP_002051.2](#)

Dilution Range WB: 1:1000, IHC-P: 1:10-50

Expiration Date 12 months from date of receipt.



Immunohistochemical analysis of formalin...