



Product Datasheet

IL1RAPL1 Antibody (orb333872)


Description	Mouse monoclonal antibody to IL1RAPL1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Purified recombinant fragment of human IL1RAPL1 (AA: 541-694) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	80kDa
Entrez	11141
Dilution Range	WB: 1/500 - 1/2000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293 using IL...



Western blot analysis of A431 (Lane1), S...



ELISA Black line: Control Antigen (100 n...