



Product Datasheet

UCP2 Antibody (orb333894)


Description	Mouse monoclonal antibody to UCP2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified recombinant fragment of human UCP2 (AA: 1-309) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	33.2kDa
Entrez	7351
Dilution Range	ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293 using UC...



ELISA Black line: Control Antigen (100 n...



Western blot analysis of human recombina...