



## Product Datasheet

### IL28A Antibody (orb333905)


<b>Description</b>	Mouse monoclonal antibody to IL28A
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, ICC
<b>Immunogen</b>	Purified recombinant fragment of human IL28A (AA: 1-200) expressed in E. Coli.
<b>Preservatives</b>	Purified antibody in PBS with 0.05% sodium azide
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	Mouse IgG1
<b>Clonality</b>	Monoclonal
<b>MW</b>	22.3kDa
<b>Entrez</b>	<b>282616</b>
<b>Dilution Range</b>	ICC: 1/200 - 1/1000; ELISA: 1/10000
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of human recombinant IL28A.



Western blot analysis of HEK293 cells using IL28A antibody.



ELISA titration curve showing O.D. vs Serial Dilutions of Antibody. Control line: Control Antigen (100 ng).