



Product Datasheet

PTPN6 Antibody (orb333907)


| | |
|----------------------------|--|
| Description | Mouse monoclonal antibody to PTPN6 |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, FC, ICC |
| Immunogen | Purified recombinant fragment of human PTPN6 (AA: 243-541) expressed in E. Coli. |
| Preservatives | Purified antibody in PBS with 0.05% sodium azide |
| Concentration | 1mg/ml |
| Storage | 4°C; -20°C for long term storage. Avoid freeze/thaw cycles. |
| Note | For research use only |
| Isotype | F |
| Clonality | Monoclonal |
| MW | 67.6kDa |
| Entrez | 5777 |
| Dilution Range | ICC: 1/200 - 1/1000; RCM: 1/200 - 1/400; ELISA: 1/10000 |
| Expiration Date | 12 months from date of receipt. |



Immunofluorescence analysis of HeLa cell...



Western blot analysis of human recombin...



Western blot analysis of HEK293 using PT...