



Product Datasheet

FOXM1 Antibody (orb333915)


Description	Mouse monoclonal antibody to FOXM1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Purified recombinant fragment of human FOXM1 (AA: 649-748) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	84.3kDa
Entrez	2305
Dilution Range	RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of human recombinant FOXM1.



Western blot analysis of HEK293 cells using FOXM1 antibody.



ELISA Black line: Control Antigen (100 nM).