



Product Datasheet

DDX39B Antibody (orb333977)


Description	Mouse monoclonal antibody to DDX39B
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Purified recombinant fragment of human DDX39B (AA: 1-250) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	49kDa
Entrez	7919
Dilution Range	WB: 1/500 - 1/2000; IHC: 1/200 - 1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of human recombinant...



Western blot analysis of HEK293 using DD...



Western blot analysis of HEK293 (Lane1),...