



Product Datasheet

RAD18 Antibody (orb334010)


Description	Mouse monoclonal antibody to RAD18
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Purified recombinant fragment of human RAD18 (AA: 332-430) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	56.2kDa
Entrez	56852
Dilution Range	RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of human recombinant RAD18.



Western blot analysis of HEK293 cells using RAD18 antibody.



ELISA Black line: Control Antigen (100 ng).