


Product Datasheet

NDRG1 antibody (Biotin) (orb334069)

Description	Goat polyclonal antibody against NDRG1 conjugated to Biotin.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	NDRG1, Biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 48kDa band observed in Human Kidney lysates (calculated MW of 42.8kDa according to NP_006087.2). See non-biotinylated parental product's datasheet for further QC data. Recommended concentration: 2-4µg/ml.
Protein Sequence	GNSAGPKSMEVSC.
Clonality	Polyclonal
MW	42.8; 35.5; 33.7; 44.7; 22.1
Entrez	10397
Dilution Range	ELISA: 1:32000, WB: 2-4 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Kidney ly...