


Product Datasheet

AIF1 antibody (Biotin) (orb334079)

Description	Goat polyclonal antibody against AIF1 conjugated to Biotin.
Species/Host	Goat
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	AIF1 / IBA1 (isoform 1 and 3), biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 17kDa band observed in Rat Brain lysates (calculated MW of 16.7kDa according to NP_001614.3). See non-biotinylated parental product's datasheet for further QC data. Recommended concentration: 0.1-0.3µg/ml.
Protein Sequence	TGPPAKKAISELP.
Clonality	Polyclonal
MW	10.5; 14.6; 16.7; 20.9; 18.8
Entrez	199
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat Brain lysat...