


Product Datasheet

GAPDH antibody (Biotin) (orb334083)

Description	Goat polyclonal antibody against GAPDH conjugated to Biotin.
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	GAPDH (C Terminus) Loading Control, biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 35kDa band observed in lysates of cel line HEK293 (calculated MW of 36.1kDa according to NP_002037.2). See non-biotinylated parental product's datasheet for further QC data. Recommended concentration: 0.3-1µg/ml.
Protein Sequence	HQVSSDFNSDT.
Clonality	Polyclonal
MW	36.1; 31.5
Entrez	2597
Dilution Range	ELISA: 1:64000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293 lysate u...