

Product Datasheet

Human Hemoglobin Subunit Theta-1 protein (orb334274)

Description	Recombinant Human Hemoglobin subunit theta-1 is produced by our E.coli expression system and the tar...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.0.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	N-His
Note	For research use only
Application notes	6His tag at the N-terminus.
Uniprot ID	P09105
Expiration Date	6 months from date of receipt.