

Product Datasheet

Human IL-23 Recetor protein (orb334285)

| | |
|--------------------------|---|
| Description | Recombinant Human Interleukin-23 Receptor is produced by our Mammalian expression system and the tar... |
| Endotoxins | Less than 0.1 ng/ug (1 EU/ug) |
| Form/Appearance | Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4. |
| Storage | Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months. |
| Tag | C-Fc |
| Note | For research use only |
| Application notes | Fc tag at the C-terminus. |
| Uniprot ID | Q5VWK5 |
| Expiration Date | 6 months from date of receipt. |