


Product Datasheet

Prion protein PrP/PRNP Antibody (orb334515)

Description	Prion protein PrP/PRNP Antibody	130KD —
Species/Host	Rabbit	70KD —
Reactivity	Human, Mouse, Rat	55KD —
Conjugation	Unconjugated	35KD —
Tested Applications	IHC, WB	25KD —
Immunogen	E.coli-derived human PRNP recombinant protein (Position: S143-S230).	15KD —
Form/Appearance	Lyophilized	WB analysis of PRNP using anti-PRNP anti...
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.	IHC analysis of PRNP using anti-PRNP ant...
Note	For research use only	
Application notes	Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.	IHC analysis of PRNP using anti-PRNP ant...
Isotype	Rabbit IgG	
Clonality	Polyclonal	
MW	27661 MW	
Uniprot ID	P04156	
Expiration Date	12 months from date of receipt.	