
Product Datasheet

SULT2B1 Antibody (orb334554)

Description

SULT2B1 Antibody

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested Applications

FC, ICC, IF, IHC, WB

Immunogen

A synthetic peptide corresponding to a sequence in the middle region of human SULT2B1 (190-218aa FDHIKGWLRMKGKDNFLFITYEELQDQLQ), different from the related mouse sequence by six amino acids, and from the related rat sequence by five amino acids.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results.
 By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users.
 . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW


41308 MW

Uniprot ID
[O00204](#)
Sensitivity


> 5000 cells

Expiration Date


12 months from date of receipt.



Flow Cytometry analysis of A549 cells us...



WB analysis of SULT2B1 using anti-SULT2B...



IF analysis of SULT2B1 using anti-SULT2B...