

Product Datasheet

Smac/Diablo Antibody (orb334593)

Description

Smac/Diablo Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ELISA, FC, ICC, IF, IHC, WB

Immunogen

E. coli-derived human Smac/Diablo recombinant protein (Position: A56-D239). Human Smac/Diablo shares 88.6% amino acid (aa) sequence identity with mouse Smac/Diablo.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality


Polyclonal

MW


27131 MW

Uniprot ID
[Q9NR28](#)
Expiration Date


12 months from date of receipt.



Flow Cytometry analysis of K562 cells us...



WB analysis of Smac/Diablo using anti-Sm...



IF analysis of Smac/Diablo using anti-Sm...