

Product Datasheet

FCGRT Antibody (orb334599)

Description	FCGRT Antibody	130KD — 100KD —
Species/Host	Rabbit	70KD —
Reactivity	Human	55KD —
Conjugation	Unconjugated	35KD —
Tested Applications	IHC, WB	25KD —
Immunogen	E. coli-derived human FCGRT recombinant protein (Position: Q56-E291). Human FCGRT shares 70.2% and 68.9% amino acid (aa) sequence identity with mouse and rat FCGRT, respectively.	15KD —
Form/Appearance	Lyophilized	WB analysis of FCGRT using anti-FCGRT an...
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.	IHC analysis of FCGRT using anti-FCGRT a...
Note	For research use only	
Application notes	Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.	
Isotype	Rabbit IgG	
Clonality	Polyclonal	
MW	39743 MW	
Uniprot ID	P55899	
Expiration Date	12 months from date of receipt.	