

Product Datasheet

ZBTB7A Antibody (orb334644)

Description	ZBTB7A Antibody	
Species/Host	Rabbit	
Reactivity	Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	IHC, WB	WB analysis of ZBTB7A using anti-ZBTB7A ...
Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of mouse ZBTB7A (125-163aa DLLERQILAADDVGDASQPDGAGPTDQRNLLRAKEYLEF), different from the related human sequence by eleven amino acids, and from the related rat sequence by one amino acid.	
Form/Appearance	Lyophilized	
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.	IHC analysis of ZBTB7A using anti-ZBTB7A...
Note	For research use only	
Application notes	Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.	IHC analysis of ZBTB7A using anti-ZBTB7A...
Isotype	Rabbit IgG	
Clonality	Polyclonal	
MW	60281 MW	
Uniprot ID	O88939	
Sensitivity	> 5000 cells	
Expiration Date	12 months from date of receipt.	