


Product Datasheet

Human RANTES protein (orb334647)

Description	Human RANTES protein
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	RANTES
Form/Appearance	20mM Tris, 150mM NaCl, 5% Glycerol, 1mM DTT pH 7.5
Concentration	0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Application notes	N-Trx-his6
MW	26 kDa
NCBI	002985
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human RANTES protei...