
Product Datasheet

Histone H2A (Phospho-S129) antibody (orb334665)

Description

Rabbit polyclonal antibody to Histone H2A

kDa A



Western blot analysis of Histone H2A (Ph...

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

WB

Immunogen

KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S129 of human Histone H2A protein. The exact sequence is proprietary.

Target

HIST1H2AG; HIST1H2AI; HIST1H2AK; HIST1H2AL; HIST1H2AM; HIST2H2AA3; HIST2H2AA4; HIST3H2A

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID
[Q8BFU2](#), [Q7L7L0](#), [P02262](#), [P22752](#), [Q6FI13](#), [Q4FZT6](#), [P0C0S8](#), [P0CC09](#), [Q6GSS7](#)
Entrez
[365877](#), [8332](#), [8330](#), [723790](#), [15267](#), [64646](#), [319162](#), [8337](#), [8329](#), [8969](#), [8336](#), [319164](#), [92815](#)
Dilution Range

WB: 1:1000-1:2000

Expiration Date

12 months from date of receipt.