

Product Datasheet

Histone H2A.X (Phospho-T120) antibody (orb334668)

Enables Biocompare.
Description

Rabbit polyclonal antibody to Histone H2A.X



Western blot analysis of Histone H2A.X (...)

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T120 of human Histone H2A.X protein. The exact sequence is proprietary.
Target	H2AFX
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P16104 , P27661
Entrez	15270 , 3014
Dilution Range	WB: 1:1000-1:2000
Expiration Date	12 months from date of receipt.