

## Product Datasheet

### Hepatitis B Virus Surface antibody (orb334812)

<b>Description</b>	Mouse monoclonal antibody to Hepatitis B Virus Surface
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Recombinant Hepatitis B Virus Surface Antigen
<b>Preservatives</b>	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.
<b>Form/Appearance</b>	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>Source</b>	Mouse
<b>Dilution Range</b>	ELISA: Use at an assay dependent dilution.
<b>Expiration Date</b>	12 months from date of receipt.