

Product Datasheet

Anti-Hepatitis B Virus E Antigen Antibody-FITC labeled (orb334859)

| | |
|----------------------------|--|
| Description | Goat polyclonal antibody against Hepatitis B Virus E conjugated to FITC |
| Species/Host | Goat |
| Reactivity | Virus |
| Conjugation | FITC |
| Tested Applications | ELISA |
| Immunogen | Recombinant Hepatitis B Virus E Antigen |
| Preservatives | Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide. |
| Form/Appearance | Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide. |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| Source | Goat |
| Dilution Range | ELISA: Use at an assay dependent dilution. |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)