

## Product Datasheet

### Hepatitis B Virus E antibody (Biotin) (orb334863)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Mouse monoclonal antibody against Hepatitis B Virus E conjugated to Biotin                       |
| <b>Species/Host</b>        | Mouse  |
| <b>Reactivity</b>          | Virus  |
| <b>Conjugation</b>         | Biotin   |
| <b>Tested Applications</b> | ELISA  |
| <b>Immunogen</b>           | Recombinant Hepatitis B Virus E Antigen  |
| <b>Preservatives</b>       | Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.             |
| <b>Form/Appearance</b>     | Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.             |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles. |
| <b>Note</b>                | For research use only  |
| <b>Clonality</b>           | Monoclonal   |
| <b>Source</b>              | Mouse  |
| <b>Dilution Range</b>      | ELISA: Use at an assay dependent dilution.   |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |