

Product Datasheet

His-tag antibody (Biotin) (orb334865)

Description	Mouse monoclonal antibody against His-tag conjugated to Biotin
Species/Host	Mouse
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	KLH-conjugated synthetic peptide HHHHHH.
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.
Form/Appearance	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.
Concentration	p>1 mg/ml
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Dilution Range	ELISA: Use at an assay dependent dilution.
Expiration Date	12 months from date of receipt.