

Product Datasheet

Hepatitis B Virus E antibody (Biotin) (orb334870)

| | |
|----------------------------|--|
| Description | Goat polyclonal antibody against Hepatitis B Virus E conjugated to Biotin |
| Species/Host | Goat |
| Reactivity | Virus |
| Conjugation | Biotin |
| Tested Applications | ELISA |
| Immunogen | Recombinant Hepatitis B Virus E Antigen |
| Preservatives | Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide. |
| Form/Appearance | Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide. |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| Source | Goat |
| Dilution Range | ELISA: Use at an assay dependent dilution. |
| Expiration Date | 12 months from date of receipt. |