

Product Datasheet

Mouse 4-1BB Ligand / TNFSF9 Protein (orb334876)

Description Mouse 4-1BB Ligand, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Arg 104 - Glu ...

Reactivity Mouse

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target 4-1BB Ligand

Preservatives Tris with Glycine, Arginine and NaCl, pH7.5

Form/Appearance Powder

Storage -20°C

Tag C-Fc

Note For research use only

Application notes This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 49.0 kDa. The protein migrates as 64-80 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence NP_033430

Purity 95%


Source Mouse 4-1BB Ligand, Fc Tag (orb334876) is expressed from human 293 cells (HEK293). It contains AA Arg 104 - Glu 309 (Accession # NP_033430).

MW 49.0 kDa

Uniprot ID [P41274](#)


NCBI [NP_033430](#)

Expiration Date 6 months from date of receipt.



SDS-PAGE analysis of Mouse 4-1BB Ligand

...



Line graph illustrates about the sensitivity...