


Product Datasheet

Human B7-H5 / Gi24 / VISTA Protein (orb334878)

Description	Human B7-H5, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Phe 33 - Ala 194 (Acc...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	B7-H5
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 44.8 kDa. The protein migrates as 60-66 kDa under reducing (R) condition (SDS-PAGE) due to Glycosylation.
Protein Sequence	AAH20568
Purity	95%
Source	Human B7-H5, Fc Tag (orb334878) is expressed from human 293 cells (HEK293). It contains AA Phe 33 - Ala 194 (Accession # AAH20568).
MW	44.8 kDa
Uniprot ID	Q9H7M9
NCBI	AAH20568
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human B7-H5 protein