


Product Datasheet

Rat PD-1 / PDCD1 Protein (orb334915)

Description	Rat PD-1, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Leu 25 - Gln 167 (Access...
Reactivity	Rat
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	PD-1
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 42.9 kDa. The protein migrates as 55-66 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_001100397.1
Purity	95%
Source	Rat PD-1, Fc Tag (orb334915) is expressed from human 293 cells (HEK293). It contains AA Leu 25 - Gln 167 (Accession # NP_001100397.1).
MW	42.9 kDa
Uniprot ID	D3ZIN8
NCBI	NP_001100397.1
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Rat PD-1 protein