



Product Datasheet

Human PD-L2 / B7-DC Protein (orb334918)

Description	Human PD-L2, Mouse IgG1 Fc Tag is expressed from human 293 cells (HEK293). It contains AA Leu 20 - P...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	PD-L2
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-mFc(mIgG1)
Note	For research use only
Application notes	This protein carries a mouse IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 48.9 kDa. The protein migrates as 60-70 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	AAI13679
Purity	95%
Source	Human PD-L2, Mouse IgG1 Fc Tag (orb334918) is expressed from human 293 cells (HEK293). It contains AA Leu 20 - Pro 219 (Accession # AAI13679).
MW	48.9 kDa
Uniprot ID	Q9BQ51
NCBI	AAI13679
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human PD-L2 protein



Line graph illustrates about the sensitivity...