


Product Datasheet

Human P-Selectin / CD62P Protein (orb334926)

Description	Human P-Selectin, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Trp 42 - Ala 771...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	P-Selectin
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 106.6 kDa. The protein migrates as 130-180 kDa under reducing (R) condition (SDS-PAGE) due to Glycosylation.
Protein Sequence	AAN06828
Purity	95%
Source	Human P-Selectin, Fc Tag (orb334926) is expressed from human 293 cells (HEK293). It contains AA Trp 42 - Ala 771 (Accession # AAN06828).
MW	106.6 kDa
Uniprot ID	P16109
NCBI	AAN06828
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human CD62P protein