


Product Datasheet

HIF1 alpha antibody (orb334997)

Description	Goat polyclonal to alpha subunit of transcription factor hypoxia-inducible factor-1 (HIF-1), which i...
Species/Host	Goat
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Purified recombinant peptide derived from within residues 470 to 590 aa of human HIF1a produced in E. coli.
Target	HIF1a
Preservatives	PBS, 20% glycerol and 0.05% sodium azide
Form/Appearance	Polyclonal antibody supplied as a 100 µl (3 mg/ml) aliquot in PBS, 20% glycerol and 0.05% sodium azide. This antibody is epitope-affinity purified from goat antiserum.
Concentration	3 mg/ml
Storage	Store at -20°C for long-term storage. Store at 2-8°C for up to one month.
Note	For research use only
Application notes	The antibody solution should be gently mixed before use.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2,000
Expiration Date	12 months from date of receipt.



Western blot analysis of staining of HT1...