



Product Datasheet

GPR55 antibody (orb335562)

Description	Rabbit polyclonal antibody to GPR55
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, ELISA: 1:20000-1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293 cells us...



Immunohistochemical staining of paraffin...