


## Product Datasheet

### E2F1 (Phospho-Ser364) antibody (orb336036)

<b>Description</b>	Rabbit polyclonal antibody to E2F1 (Phospho-Ser364)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	WB: 1:1000-3000
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of EGF treated 293...