



## Product Datasheet

### FGFR3 (Phospho-Tyr577) antibody (orb336101)


<b>Description</b>	Rabbit polyclonal antibody to FGFR3 (Phospho-Tyr577)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	WB: 1:1000-3000
<b>Expiration Date</b>	12 months from date of receipt.



Immunohistochemical staining of Mouse testis



Immunohistochemical staining of Rat testis



Line graph illustrates about the Ag-Ab r...