



Product Datasheet

LTK/ALK (Phospho-Tyr672/Tyr1278) antibody (orb336105)

Description	Rabbit polyclonal antibody to LTK/ALK (Phospho-Tyr672/Tyr1278)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:1000-3000, ELISA: 1:20000-40000
Expiration Date	12 months from date of receipt.



Western blot analysis of H2O2 treated K5...



Western blot analysis of H2O2 treated K5...