

Product Datasheet

Cdk1/2 (Phospho-T14) antibody (orb336535)

Description Rabbit polyclonal antibody to Cdk1/2 (Phospho-T14)

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IHC-P, WB


Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Clonality Polyclonal

Dilution Range WB: 1:1000-3000, IHC-P: 1:200

Expiration Date 12 months from date of receipt.



Western blot analysis of 293 whole cell ...